



NEW GENERATION DEVICES

The Body Knows the Difference®

Flat Osteotomy Blades

Some Indications for Use Include:

- THR
- Mandibulectomy
- Rhinotomy
- PAUL
- Corrective osteotomy
- SHO
- Hemipelvectomy
- Femoral shaft transaction during amputation
- Arthrodesis cuts
- Sagittal osteotomy



10mm



20mm



30mm

The NGD flat osteotomy blades have been designed for the RIGID Job Max® multitool.

These blades can be used with other oscillating multitool brands with the addition of a universal adapter.

The current list of compatible multitools includes; Chicago, Craftsman, Dremel, Electric, Fein 250 Series, Mastercraft, Masterforce, Performax, Rockwell and Tool Shop.

***Reach for a better standard of care
with confidence***

***Contact us today to place your
order, or learn more.***

sales@ngdvet.com

1 800 797 8587